

Carboxyl QDs conjugation protocol

Material provided	Storage
Quantum dots nanoparticles with carboxyl (-COOH) functionalisation, 1 mL at a concentration of 2 μM.	- Store at 2–8 °C. - Do not freeze! - Protect from light.

Material NOT provided
- Antibody solution at 1 mg/mL in a suitable buffer (minimum antibody concentration: 0.5 mg/mL), such as Dubecco's PBS (DPBS) and HEPES. <i>Amine-containing buffers such as Tris are <u>not</u> suitable for this reaction!</i> 1-Ethyl-3-(3-dimethylaminopropyl)carbodiimide (EDC). - N-hydroxysuccinimide (NHS) or N-hydroxysulfosuccinimide (sulfo-NHS). - Glycine (for quenching). - Centrifugal filter with a MWCO of 150 kDa, such the Amicon® Ultra 0.5 mL (Merck Millipore).

Procedure to conjugate 100 μ g of antibody to QD-COOH

Step 1 – QDs activation

1. In a 1.5 mL tube, add 140 μ L of carboxyl QDs to 140 μ L of DPBS.
2. Freshly prepare a 10 mg/mL solution of EDC and a 12 mg/mL solution of NHS, both in ultrapure water.
3. Immediately add 10 μ L of the EDC solution to the QDs solution, followed quickly by 10 μ L of the NHS solution.
4. Mix the QDs at room temperature for 20 minutes.

Step 2 – QDs + antibody reaction

1. Purify the activated QDs from excess EDC and NHS via spin filtration, repeating the filtration three times and resuspending the nanoparticles in DPBS after each step.
NOTE: it is important to complete this step as quickly as possible, as activated QDs will hydrolyse and lose their reactivity.
2. Collect the purified QDs, preferably in a 200 μ L tube. Adjust the final volume to 100 μ L using DPBS.

3. Add 100 μ L of the antibody solution to the QDs, and mix for at least 3 hours at 4 °C (overnight incubation is also acceptable).
4. Quench any unreacted NHS groups by adding 15 μ L of 50 mM glycine.
5. Store the QD–antibody conjugate at 2–8 °C, protected from light. Do not freeze! The conjugate should remain stable for at least 1 month.

NOTE: To conjugate a smaller amount of protein, take an aliquot of the carboxyl-functionalised QDs and place it in a new tube. Maintain the same molar ratios between EDC, NHS, antibody, buffer, and QDs.

NOTE: Removal of excess antibody can be achieved via size exclusion chromatography or centrifugal filtration with a MWCO > 150 kDa.