

# Azide QDs conjugation protocol

Material provided	Storage
<ul> <li>Quantum dots nanoparticles with azide (-N<sub>3</sub>) functionalisation, 500 µL at a concentration of 2 uM.</li> <li>Borate Buffered Saline (BBS) buffer, 10 mM, pH 8.4, 8 g/L NaCl</li> <li>Linker solution (0.1 mM)</li> </ul>	<ul> <li>Store at 2–8 °C.</li> <li>Do not freeze!</li> <li>Protect from light.</li> </ul>

### Material NOT provided

- Antibody solution at 1 mg/mL in a suitable buffer (minimum antibody concentration: 0.5 mg/mL), such as BBS, DPBS, and HEPES. Amine-containing buffers such as Tris are not suitable for this reaction!
- Dubecco's PBS (DPBS).
- Glycine (for quenching).
- Centrifugal filter with a maximum MWCO of 50 kDa, such the Amicon® Ultra 0.5 mL (Merck Millipore).

## Procedure to conjugate 100 $\mu$ g of antibody to QD-N<sub>3</sub>

#### Day 1 – Linker + antibody reaction

- 1. In a 1.5 mL tube, add 100  $\mu L$  of antibody (1 mg/ml) to 500  $\mu L$  of BBS.
- Immediately add 12 µL of the linker solution to the tube containing the antibody solution and mix for 18 hours at 4 °C.
   (These quantities give the preferred molar ratio of linker to antibody of 2:1.)

#### Day 2 – Linker-antibody + QDs reaction

- 1. Quench the linker reaction with  $15 \,\mu\text{L}$  of 50 mM glycine.
- 2. Purify the antibody from excess linker via spin filtration three times, using DPBS to resuspend the protein after each step.
- 3. Collect the concentrated protein, preferably in a 200  $\mu L$  tube. Add enough DPBS to adjust the final volume to 100 uL.
- 4. Add 140 µL of QDs to the antibody-linker solution, and mix for 3 hours at 4 °C.
- 5. After 3 hours, the solution is ready to use.



6. Store the QD-antibody conjugate at 2–8 °C, protected from light. Do not freeze! The conjugate should remain stable for at least 1 month.

<u>NOTE</u>: To conjugate a smaller amount of protein, take an aliquot of the azidefunctionalised QDs and place it in a new tube. Maintain the same molar ratios between the linker solution, antibody, buffer, and QDs. If the scale of the reaction requires less than 1  $\mu$ L of linker, dilute the linker 1:10 in DMSO (final concentration 0.01 mM), and immediately use the diluted solution for the antibody modification.